

PHILIPS formerly ADAC

CAMERA MODEL	CONFIGURATION	CAMERA UFOV <small>USEFUL FIELD OF VIEW</small>	SOURCE ACTIVE DIMENSIONS	PERFLEXION PRODUCT CODE	FEATHER LITE MODEL	C-THRU MODEL	CARDINAL	MODALITY
411	Circular	11in	14in	N/A	3700.AD.005M.N	N/A	N/A	SPECT
412	Circular	12in	14in	N/A	3700.AD.005M.N	N/A	N/A	SPECT
415	Circular	15in	18.5in	PF18C-057-5M	3706.AD.005M.N	CTRF10051	SFL18C-057-5m	SPECT
ARC 3000	Circular	15in	18.5in	PF18C-057-10M	3706.AD.010M.N	CTRF10050	SFL18C-057-10M	SPECT
Argus	Rectangular	20in x 15in	23.9in x 16.4in	PF24R-057-10M	3709.AD.010M.N	CTRF10000	SFL24R-057-10M	SPECT
Brightview XCT	Rectangular	20.25in x 16in	24.48in x 17.39in	PF24RW-057-10M	N/A	CTRF10021	N/A	SPECT
Cardio 60	Rectangular	20in x 15in	16in x 24in	PF24R-057-10M	3709.AD.010M.N	CTRF10000	SFL24R-057-10M	SPECT
Cardio MD	Rectangular	9.2in x 15.4in	16in x 10in	PF16R-057-10M	3727.AD.010M.N	CTRF10017	N/A	SPECT
Cirrus	Circular	15in	18.5in	PF18C-057-10M	3706.AD.010M.N	CTRF10050	SFL18C-057-10M	SPECT
Forte	Rectangular	20in x 15in	16in x 24in	PF24R-057-10M	3709.AD.010M.N	CTRF10000	SFL24R-057-10M	SPECT
Genesys	Rectangular	20in x 15in	23.9in x 16.4in	PF24R-057-10M	3709.AD.010M.N	CTRF10000	SFL24R-057-10M	SPECT
Precedence	Rectangular	20in x 15in	16in x 24in	PF24R-057-10M	3709.AD.010M.N	CTRF10000	SFL24R-057-10M	SPECT
Solus, Cardial, or Vertex	Rectangular	20in x 15in	16in x 24in	PF24R-057-10M	3709.AD.010M.N	CTRF10000	SFL24R-057-10M	SPECT
Skylight	Rectangular	20in x 14in	16in x 24in	PF24R-057-10M	3709.AD.010M.N	CTRF10000	SFL24R-057-10M	SPECT

PHILIPS formerly Marconi/Picker

Axis/ Irix	Rectangular	21in x 15.5in	16in x 24in	PF24R-057-10M	3709.AD.010M.N	CTRF10000	SFL24R-057-10M	SPECT
Dynamo	Circular	10in	14in	N/A	3700.AD.010M.N	N/A	N/A	SPECT
Prism 1000	Rectangular	20in x 15in	16in x 24in	PF24R-057-10M	3709.AD.010M.N	CTRF10000	SFL24R-057-10M	SPECT
Prism 2000	Rectangular	20in x 15in	16in x 24in	PF24R-057-10M	3709.AD.010M.N	CTRF10000	SFL24R-057-10M	SPECT
Prism 3000	Rectangular	15.7in x 9.4in	16in x 10in	PF16R-057-10M	3727.AD.010M.N	CTRF10017	N/A	SPECT
SX300	Square	14in x 14in	Various	N/A	MED3725	CTRF10000	SFL14R-057-10M	SPECT

Sources for Philips Scanners/Cameras

CAMERA MODEL	PRODUCT CODE	Qty	DESCRIPTION	NUCLIDE	ACTIVITY	CALIBRATION	MODALITY
Cardio 60, Vertex 60, Forte & Skylight with Vantage Pro	PHL-Forte	1	Source Set Consists of 2 ea NES8412 3 ea 078410 (Lead Pigs)	Gd-153 N/A	250mCi(9.25GBq) N/A	Nominal N/A	SPECT
Allegro	GF-0227	1	Calibration Disc	Na-22	100uCi(3.7MBq)	Nominal	PET
Brightview XCT	MMS07-153-100U	6	Multi Modal QC Disc	Gd-153	100uCi(3.7MBq)	Nominal	SPECT
CPET	GF-0227	1	Calibration Disc	Na-22	100uCi(3.7MBq)	Nominal	PET
Cardio MD	PHL-Cardio MD	1	Source Set Consists of 2 ea NES8497	Gd-153	240mCi(8.88GBq)	Nominal	SPECT
Gemini TF / LXL	HEGL-0136	1	QC Rod	Na-22	100uCi(3.7MBq)	Nominal	PET
	MMS05-022-10U	6	Multi Modal QC Disc	Na-22	10uCi(370kBq)	Nominal	PET
Gemini GXL	GF-0227	1	Calibration Disc	Na-22	100uCi(3.7MBq)	Nominal	PET
	MMS05-022-10U	6	Multi Modal QC Disc	Na-22	10uCi(370kBq)	Nominal	PET
Ingenuity PET/MR	RFQ1021-022-20U	6	PET / MR Disc	Na-22	20uCi(740kBq)	Nominal	PET
	HEGL-0136	1	QC Rod	Na-22	100uCi(3.7MBq)	Nominal	PET
Ingenuity TF	HEGL-0136	1	QC Rod	Na-22	100uCi(3.7MBq)	Nominal	PET
	MMS05-022-10U	6	Multi Modal QC Disc	Na-22	10uCi(370kBq)	Nominal	PET
Mosaic Animal PET	GF-0227	1	Calibration Disc	Na-22	100uCi(3.7MBq)	Nominal	PET